


How To Hook Up Logger Pro Motion Sensors



- LABPRO SHOULD BEEP WHEN POWER IS PLUGGED IN.
- WHEN ALL THE EQUIPMENT IS HOOKED UP, OPEN "LOGGER PRO"
 - FIND IT UNDER "ALL PROGRAMS" IN THE "VERNIER" FOLDER
- TO RECORD MOTION, HIT THE COLLECT BUTTON (IT WILL BE GREEN)

POSITION

Time

VELOCITY

Time

ACCELERATION

Time